



Vein Finder

VS30

NEW
Model

VEIN FINDER

Visualization Modes

Brightness

Adjustable brightness allows clinicians customizing the image brightness base on the environment light.

Colors

Switchable colors can fit different needs from the nurses.

Size

Different image size can help medical staff focus on one veins, which is also useful for child.

Inverse

Clinicians can highlight the veins or highlight the surrounding tissue.

VS30 Parameter

Light Type	Near-infrared light
Best Projection Distance	200±20mm
Dimension	240x60x60mm(LxWxH)
Battery Working Time	4 hours
Operating Mode	Handheld or suspension

VS30 Optional Stand



Main Unit ▶



VS30 Application:

vein finder can digitally displays a map of the vasculature on the surface of the skin in real time, aiding healthcare professionals to verify vein patency and avoid valves or bifurcations. Vein finder is a proven modality that can improve patient care and save valuable time for both clinicians and patients.


Adult



Forehead(Child)



Obesity Patient(hand)



Head



Dark Skin(hand)

