

CE

WYKL5

Shadowless Operating Light




User-friendly and Best Value for Money



CAD/CAM designed multi-angle reflectors, eliminate shadows from the surgeon's heads and hands



Switch on/off switch in supporting tube



Light focusing adjustment knob



Brightness adjustment in the base

WYKL5

Shadowless Operating Light

Premium Quality, Surprising Price

WYKL5

Shadowless Operating Light

General Features

High-intensity of over 80,000 lux at 1 meter;

4 000K color temperature offers natural and true coloring;

Uses five halogen bulbs – each rated with a 2000-hour bulb life; bulbs wired in parallel for continuous operation;

Extremely smooth focus mechanism with centre handle;

Cool operation for long hours of use;

Two silence brakes for safety use;

Steady and smooth movement of the light head;

Durable and cost-efficient.

Technical Specifications

Illumination: $\geq 80,000$ lux, adjustable

Color temperature: 4000 ± 500 K

Color rendering index: 85%

Light field size: 150 ± 35 mm

Illumination depth: 450mm

Qty of light bulbs: 5 pcs

Bulb power: 24V/25W


Diameter of lamp head: 500mm

Temperature rise at surgeon's head: Less than 10°C

Optional backup battery: UPS, 24ah, last around 3 hours once fully charged. (illumination fixed with UPS)



Ref picture for WYK5 with backup battery



All the measurement is millimetre.