

# TD4Z

## Portable Low-Speed Centrifuge



### 1 Technical performance

- It is popular choice in any clinic and laboratory.
- Microcomputer control, brushless motor.
- It is programmable and digital display.
- It has speed/RCF switch easily.
- There is lid safety device, speed safety device, and error detection.
- Silastic sealed ring comfort to GMP.US FDA certificated.

### 3 Rotor



No.1

Capacity: 12x7ml  
Max Speed: 4000r/min  
Max RCF: 1980xg



No.2

Capacity: 8x15ml  
Max Speed: 4000r/min  
Max RCF: 1980xg



No.3

Capacity: 8x15ml  
Max Speed: 4000r/min  
Max RCF: 1980xg



No.3

Capacity: 8x20ml  
Max Speed: 4000r/min  
Max RCF: 1980xg



No.3

Capacity: 12x10ml  
Max Speed: 4000r/min  
Max RCF: 1980xg

### 2 Specification

Model	TD4Y
Max speed	4000r/min
Max RCF	1980xg
Max capacity	8x20ml
Speed accuracy	± 30r/min
Time setting	1min-99min
Noise	<50dB(A)
Power supply	AC110V - AC220V
Consumption	80W
Dimension(L×W×H)	310mm×260mm×220mm
Packing size(L×W×H)	370mm×320mm×280mm
Net weight (No rotor)	7kg