



LW W12

Microplate Washer

- 12-way wash head interchangeable
- Programmable up to 64 wash protocols
- Minimal residual volume
- Low well-to-well variation
- Automatic self-diagnostic program
- Automatic liquid monitoring system
- Touch screen , LCD display

LANDWIND
MEDICAL

• Specifications



- **WASH HEAD**
12 or 8 heads
- **WASH MODE**
Row and plate
- **WASH ROWS**
1-12 rows
- **PRIME VOLUME**
50-700 μ l
- **MOVING CYCLE**
1-9 cycles
- **CLEANING TIME**
1-255s
- **CLEAR BOTTOM**
Yes/No
- **WASH BOTTOM**
Yes/No
- **BOTTLE CAPACITY**
>2L
- **WASH PROGRAMS**
Up to 64
- **TYPE OF PLATE**
Up to 10
- **RESIDUAL VOLUME PER WELL**
<2 μ l for V/ U bottom plates
<3 μ l for flat bottom plates
- **DISPLAY**
Backlight LCD touch screen
- **KEYPAD**
Touch screen

- **OPERATING ENVIRONMENT**

| | |
|---------------------|---------------------------------|
| Input power | 80W |
| Power supply | AC220V/110V, 50/60Hz |
| Ambient temperature | 10 $^{\circ}$ C-30 $^{\circ}$ C |

| | |
|-------------------|--|
| Relative humidity | 20%-85% |
| Weight | 7.5Kg (Bottles not included) |
| Dimension | 380mm(H) \times 350mm(W) \times 190mm(D) |

Shenzhen Landwind Industry Co., Ltd

Headquarter:
Landwind Science & Technology Park, Tangtou Avenue,
Shiyan Town, Bao'an District, Shenzhen China 518108
Tel: +86-755-83933788 Fax: +86-755-27353240
Email: Sales@landwindmedical.com
Web: www.landwindmedical.com

DISTRIBUTOR:

Specifications and appearance are subject to change without prior notice.

P/N: 70150285

2013.06