



Examination lamp

Product Features

- *Mini lamp holder with diameter of 75mm;
- *Illumination using Osram imported from Germany, with stable light and long working time;
- *Flexible arm can freely position,easy to move;
- *Even bright white light.



Technical parameters

| | |
|--------------------|------------|
| Model | JSF-JC02 |
| Bulb | LED |
| Supply Voltage | 220V |
| Power | 3W |
| Color temperature | 4500K |
| Illumination(50cm) | 16000Lux |
| Spot size adjust | × |
| Brightness adjust | × |
| Height adjust | √ |
| Lifetime | 20000h |
| Packing size | 93*48*10cm |
| N.W/G.W | 4KG/5KG |